



SILUET
DOUBLE-UP

Find your hidden facial
line
with siluet

Diode Laser

PRINCIPLE


SILUET DOUBLE UP has two wavelengths ,980nm & 1470nm with each fiber treatment.

It is very innovative device combined double fiber each very useful wavelengths.

980nm would be very helpful to proceed a lifting effect because it has very high water absorption and very low absorption of fatty tissue.

1470nm has very high absorption for water & fatty tissue.


It would be definitely the best solution for lipolysis and precise facial & body contouring.



[980nm Laser]

Lifting Mode

Removal of small volumes of fat with
Concurrent subdermal tissue contraction



[1,470nm laser]

Subcutaneous fat layer Lipolysis

FEATURES

Minimal invasive lipolysis laser with safe procedure without scar

Optimal facial & body contouring and skin lifting at once

No downtime after treatment

Available for make up, face washing and shower within treatment day

Fast treatment(within 30min) & Long lasting result

Selective 2 fiber options(400nm, 600nm) according to target width

1kHz duration allows minimal risk for thermal damage

INDICATION

Face

Forehead Line

Frown Line

Eye Bag

Cheek born

Nasolabial Fold

Vertical Lip Lines

Cheek bone

Jowl Line

Body

Forearm

Breast Lifting

Upper abdomen

Flank

Lower abdomen

Hip Lifting

Thigh

Stretch mark

Calf

SYSTEM&ACC



Handpiece



Fiber



SPECIFICATIONS

Irradiation Type	Diode Laser
Wavelength	980nm, 1470nm
Maximum Energy	30W, 15W
Mode	CW mode, Pulse mode
Repetition Rate	1~1000Hz
Delivery	400 μm , 600 μm Optical Fiber