


**3-MAX** coolShaping™

## Cryo Sculptor Handpiece

Removal of fat cells with Cryolipolysis, Vacuum,  
and LED technology

1~5 step vacuum level adjustment, various pulse  
mode available.




## S-Cavity Handpiece

Powerful cavitation through 36KHz ultrasonic energy  
effectively reduces fat cells.

## S-Polar Handpiece

Combination of Multi-polar RF and LED helps to improve elasticity  
3 Size of handpieces allow to treat various areas from face to  
whole body.



## S-Polar Handpiece



Large



Medium



Small

## S-Sound Handpiece

1MHz Ultrasound Frequency generates heat through vibration, induces muscle relaxation, making it effective in improving cellulite.